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Other material that are not provided by this kit:

- Samplers and 1/5mL microtubes

Sample: Before using the kit, you need homogenize the Plant samples. Harvested tissues should be used freshly or stored at very low temperature as quickly as possible. Grinding in mortar and pestle under liquid nitrogen is a good method for homogenizing the sample, but alternative methods, such as a homogenizer or a bead-beater, can be employed in case by case for efficient disruption. Shaking or vortex during incubation for lysis may greatly accelerate the efficiency of lysis, resulting in reduced time for complete lysis. Note that the freshness and the particle size of disrupted sample is the key for good result and that the frozen sample should be kept on ice until use. To homogenize the plant samples, add the harvested samples (100-200 mg) to the sterile tubes and add 500 μ L Homogenize buffer (CN: KPG-HBp) and then homogenize the samples and then follow the next protocol.

Protocole:

- Add 500 microliters of lysis buffer A to a DNase/RNase free microtube.
- Add the appropriate amount of the sample (for example, 200 microliters of homogenized plant tissue) to the microtube and mix vigorously for 15 seconds on the vortex and incubate at room temperature for 1 minutes. Do not increase the incubation time, because it is harmful for DNA.
- Add 500 microliters of lysis buffer B to the microtube. And vortex for completely.
- Add 100 μ L of Lysis Buffer C to the microtube and mix vigorously for 5 seconds.
- Add 500 microliters of cold precipitating buffer to the microtube and mix gently by hand. Pay attention, you can mix lysis solution C and precipitation buffer together and then add it to the microtube.
- Add 750 μ L of the contents of the microtube to the DNA extraction columns and centrifuge at 12.000 RPM for 1 minute. Make sure that according to the kit method, very good amounts of DNA are separated and you should not add extra sample volume to the column. But if your sample is small, you should use another column for the remainder (an additional column is not provided in the kit and you have to get it separately from Karmania Pars Gene Company).
- Transfer the column to a new collection tube and add 500 μ L of washing buffer and then centrifuge at 8.000 RPM for 1 minute.
- Repeat the washing process with the washing buffer one more time.
- After finishing the wash without changing the collection tube, centrifuge the column at 12/000 RPM for 1 minute to drain the remaining buffer from the filters.
- Transfer the column to a sterile RNase/DNase free 1/5 microtube and add 100 microliters of DNase/RNase free distilled water and incubate for 2 minutes at room temperature. Make sure that the temperature of the distilled water is around 55 °C for better DNA separation.
- Next, centrifuge the columns at 12/000 RPM for 1 minute. The resulting solution contains DNA. To increase the amount of extracted DNA, we suggest transferring the extracted DNA sample to the column and centrifuging again at 12.000 RPM for 1 minute.

<http://kpgene.ir>

09132926113

@karmaniaparsgene1

info@kpgene.ir

karmaniaparsgene

Troubleshooting:

- The current kit is highly quality controlled and the probability of encountering a problem is very low. However, in some cases, you may encounter problems using this kit. Refer to the following table to solve these issues:

Problem	Reasons	Solving
Low Yield	Low cell count in the sample	✓ Use more samples
		✓ Add 30 μ L DNase/RNase free water
	Using an expired kit	✓ Although the components of the kits has more than 1 year stability, lysis buffer A and B had 6 months stability.
	Using the material contaminated with DNase	✓ Use the DNase/RNase free materials, You can prepare them from Karmania Pars Gene Company
	Using low amounts of lysis buffer	✓ Please follow the protocole
	Using column from other source	✓ Please use the column from Karmania Pars Gene Company
DNA is not appropriate for PCR	Presence of PCR inhibitor	✓ Do not change the washing step and centrifuge conditions

Safety:

The lysing solution in the kit has a very alkaline and acidic pH. On the other hand, the substances used in these solutions may in some cases be stimulating for the nervous system. Therefore, when working with these solutions, be sure to use gloves and in case of contact with the skin, wash with large amounts of water. In case of contact with the eyes, wash with large amounts of water and if there is no improvement, refer to medical centers.

KPG
Karmania Pars Gene