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Other materials needed: Samplers, Centrifuge, 1.5 mL tubes

Sample: This kit is designed to extract DNA from HPV positive samples. Using the DNA extraction kit of Karmania Pars Gene Company, you can extract large amounts of HPV DNA from all samples, even tissue. With the difference that if tissue is used, first the tissue sample must be completely homogenized in the homogenizing buffer of Karmania Pars Gene Company (CN: KPG-HBt) (the desired buffer can be prepared separately). **There is no need for this buffer when extracting from the transfer medium.** It should be noted that the current kit is able to isolate genomic DNA from blood, cell culture and other viruses. To use samples other than HPV medium, use the following values:

1. Blood, serum and plasma: 200 microliters.
2. Cell culture: up to 1×10^7 cells.
3. Homogenized tissue: 30 to 50 mg of homogenized tissue in a volume of 500 microliters of homogenizing buffer and finally take 200 microliters.

Protocol:

To extract HPV-DNA, the complete transfer medium (LBC) should be mixed first and 1.5 ml of it should be added to a 1.5 ml DNase/RNase free microtube and then centrifuged at 13.000 rpm for 3 minutes. Drain the supernatant completely and apply the following to the remaining sediment:

- Add 250 μ L of lysis buffer A to the microtube containing the sediment and mix vigorously for 5 seconds on a vortex and incubate at room temperature for 5 minutes.
- Add 250 microliters of lysis buffer B along with 10 microliters of the carrier (Carrier) to the microtube and mix vigorously for 5 seconds on the vortex and incubate at room temperature for 5 minutes.
- Add 100 μ L of lysis buffer C to the microtube and vortex vigorously for 5 seconds and incubate at room temperature for 5 minutes.
- Add 300 microliters of cold sedimentation buffer to the microtube and mix gently with hand movement and incubate at -20°C (freezer) for 5 minutes.
- Add the contents of the microtube to the DNA extraction columns and incubate for two minutes at room temperature and then centrifuge at 12.000 RPM for 30 seconds.
- Transfer the column to a new collection tube and add 500 microliters of washing buffer and then centrifuge at 12/000 RPM for 30 seconds.
- Transfer the column to a new collection tube and add 40 microliters of DNase/RNase free distilled water and incubate for 2 minutes at room temperature. Make sure that the temperature of the distilled water is around 60°C so that HPV-DNA purification can be done better.

Next, centrifuge the columns at 12/000 RPM for 1 minute. The obtained solution contains HPV-DNA and it is suggested to use this purified HPV-DNA to perform PCR in the amount of 5 μ l in the final volume of 20 μ l PCR.

Troubleshooting:

The present kit is highly quality controlled and the probability of encountering a problem is very low. However, in some cases, you may encounter problems using this kit. Refer to the following table to solve these issues:

Troubleshoot	Reason	Solving
Low yeild	Low sample	Use more cell size
	Expired kit	Use the kit with valid ET
	Existence of DNase	Use DNase free microtubes and sampler tips
	Using low volume of lysis buffer	Purify according to the kit protocol
PCR false negative	Existense of inhibitors	Washing the column one more time

Necessary precautions for working with the kit and its contents

➤ The lysing solution in the kit has a very alkaline and acidic pH. On the other hand, the substances used in these solutions may in some cases be stimulating for the nervous system. Therefore, when working with these solutions, be sure to use gloves and in case of contact with the skin, wash with large amounts of water. In case of contact with the eyes, wash with large amounts of water and if there is no improvement, refer to medical centers.

The rest of the materials in the kit are safe and do not cause damage to human tissue.

Contact us:

<http://karmaniaparsgene.com>

karmaniaparsgene@gmail.com

09132926113 +983434208024-5