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DNA extraction kit from blood and tissue

By sedimentation method (CN: KPG-DNKtbs) fifty tests

DNA extraction steps

Kit Contents:

- Lyse Buffer DNK (CN: KPG-DNK); 2 vials of 25 ml
- Lysis buffer DNC (CN: KPG-DNC); 2 vial of 10 ml
- Tissue lysis buffer (CN: KPG-LBt); 1 vial of 25 ml
- Precipitating buffer (CN: KPG-PS); 1 vial of 30 ml
- Washing buffer (KPG-WBp); 2 vials of 25 ml
- DNase free distilled water (CN: KPG-DW); 1 vial of 5 ml

Other required items:

Variable 1000 and 100 microliter sampler

Sample: This kit is designed to extract DNA from blood and tissue. Using this kit, you will be able to purify serum, plasma, VTM cell-free transfer medium, gram-negative bacteria and cell cultures, and DNA homogenized tissue samples in addition to blood. For this, use the following volumes:

1. Blood/plasma/VTM cell-free transfer medium: 200 μ l
2. Cell culture medium: at the rate of 1×10^7 cells
3. Thoroughly homogenize 30 to 50 mg of tissue in a volume of 400 μ l of tissue lysis buffer and finally use 200 μ l of the homogenized sample for DNA extraction. Be careful that the lysis buffer is harmful to DNA, so the sample should not remain in the buffer for more than 5 minutes and the next step should be started immediately.

- Add 500 microliters of lysis buffer DNK to a DNase/RNase free microtube.

- Add the appropriate amount of the sample (for example, 250 microliters of blood or 200 microliters of tissue lysate solution) to the microtube and mix vigorously for 15 seconds on the vortex and incubate at room temperature for 10 minutes. do. Vortex the samples every 5 minutes.

- Add 200 μ L of Lysis Buffer DNC to the microtube and mix vigorously for 5 seconds.

- Centrifuge the microtube at 12000-14000 RPM for 5 minutes and transfer the supernatant to another DNase/RNase free microtube.

- Add 300 microliters of cold precipitating buffer to the new microtube and mix gently by hand.

- Centrifuge the microtube at 12000-14000 RPM for 5 minutes and drain the supernatant completely.

- Add 500 microliters of washing buffer to the microtube and mix gently by hand.

- Centrifuge the microtube at 10000 RPM for 5 minutes and drain the supernatant completely and dry the remaining washing solution in the microtube by inverting it on a paper towel.

- Add 50 microliters of DNase/RNase free distilled water and incubate for 2 minutes at room temperature. Make sure that the temperature of the distilled water is around 55 degrees Celsius for better DNA separation.

Troubleshooting:

The present kit is highly quality controlled and the probability of encountering a problem is very low. However, in some cases, you may encounter problems using this kit. Refer to the following table to solve these issues:

Troubleshoot	Reason	Solving
Low yeild	Low sample	Use more cell size
	Expired kit	Use the kit with valid ET
	Existence of DNase	Use DNase free microtubes and sampler tips
	Using low volume of lysis buffer	Purify according to the kit protocol
PCR false negative	Existence of inhibitors	Washing the column one more time

Necessary precautions for working with the kit and its contents

□ The lysing solution in the kit has a very alkaline and acidic pH. On the other hand, the substances used in these solutions may in some cases be stimulating for the nervous system. Therefore, be sure to use gloves when working with these solutions and wash with large amounts of water in case of contact with the skin. In case of contact with eyes, wash with large amounts of water and if there is no improvement, refer to medical centers. The rest of the materials in the kit are safe and do not cause damage to human tissue

Contact us:

<http://karmaniaparsgene.ir>

karmaniaparsgene@gmail.com

09132926113 +983434208024-5