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Sample:

This kit is designed to extract RNA from the Sara-COV-2 virus. With this kit you will be able to discharge serum, plasma, VTM-cell-free transmission, cell patrol and hemogenized tissue RNA Sara-COV-2.

Use the following volumes for this:

1. Blood/ Plasma/ Non -Cellular VTTM: 200 MW
2. Cell Medical Media: 107 × 1 cell
3. Homogenized tissue: 30 to 50 mg hemogenized tissue in 500 microit buffer volume and finally 200 microlate harvest

Material needs

Real-Time PCR machine	Sampler tips
Centrifuge	Samplers

Protocol:

1. Add 400 microlitrs from the Liz buffer to a DNase/RNase free microtube.
2. Apply the sample (for example 200 microliter from the VTM environment) to the microtube and mix heavily on the vortex for 15 seconds and incubate for 6 minutes at room temperature.
3. Add 10 microliters from the carrier (CAR: CN: KPG-CR) to the microtube and mix severly for 5 seconds. Note that due to the nature of the carrier, which is made of temperature -resistant carbohydrate, you can keep the carrier for 1 month at ambient temperature under the hood and with respect to the conditions for microbial contamination. Keep this substance at -20 ° C if you need to maintain a carrier for more than a month.
4. Add 200 microliters of Percipitation buffer to the microtube and mix heavily by vortex for 5 seconds.
5. Transfer the microtube contents to the RNA extraction columns and centrifuge at 8.000 RPM for 1 minute.
6. Transfer the column to the new collection tube and add 500 microliters of washing buffer and then centrifuge in 8.000 RPM for 1 minute.
7. Repeat the centrifugal operation, but this time in 12.000 RPM and centrifuge for 1 minute.
8. Transfer the column to the new tube collection and add 70 microliter of DNase/RNase free water and incubate for 2 minutes at room temperature. Make sure that the distilled water temperature is about 60 degrees Celsius to better separate the SARS-COV-2-RNA.
9. Continue the columns in 12.000 RPM for 1 minute for 1 minute. The obtained solution contains Sars-COV-2-RNA and can be used for RT-PCR or stored at -20 ° C for six months and -70 ° C for more than six months.

Safety while using the kit

The solutions used in the kit have oxidizing and acidic properties. Avoid direct contact with skin and eyes. In case of contact with the mentioned tissues, wash with plenty of water and go to the nearest medical center.

Manufacture Address

Refsanjan, 20 km from Rafsanjan road to Kerman, Rafsanjan special economic zone, food and drug area
09132926113
09135025983
Karmaniaparsgene@gmail.com
www.karmaniaparsgen.ir

Other DNA extraction kits and RNA produced by Karmania Pars Gene:

All the following DNA extraction kits are available in all three columns, sedimentary and magnet nanoparticles:

1. DNA extraction kits from blood and tissue by column
2. DNA extraction kits from gram-positive and negative bacteria
3. DNA extraction kit from mycobacterium
4. DNA extraction kit from Fungi
5. DNA extraction kit from virus
6. DNA extraction kit from HPV virus
7. DNA extraction kits from plant tissue

RNA extraction kits:

All the following RNA extraction kits are available in both column and sedimentary methods:

1. RNA extraction kits from blood and tissue
2. RNA extraction kits from gram-positive and negative bacteria
3. RNA extraction kit from mycobacterium
4. RNA extraction kit from Fungi
5. The RNA extraction kit from the virus
6. RNA extraction kits from plant tissue

One of the most efficiently discharged kits: **Plasmid extraction Kit**