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karmaniaparsgene

In the name of God

Magnetic bead RNA extraction kit from blood and tissues

Kit components:

- Lyse A buffer (CN: KPG-KLSA): 1 vial 80 mL
- Lyse B buffer (CN: KPG-KLSB): 1 vial 10 mL
- Tissue lysis buffer (CN: KPG-LBt): 1 vial 25 mL
- Washing buffer (CN: KPG-WBp): 1 vial 80 mL
- RNase/DNase free water (CN: KPG-DW): 1 vial 30 mL
- Magnetic Nanoparticle (CN: KPG-Mabead): 200 mg
- 1/5 Microtubes (KPG-1/5m) 100 preps

Other material that are not provided by this kit:

- 1000 and 100 microliter samplers

Sample: This kit is designed to extract RNA from blood and tissue using magnetic nanoparticle technology. Using this kit, you will be able to purify RNA from serum, plasma, VTM cell-free transfer medium, cell culture and homogenized tissue samples in addition to blood. For this, use the following volumes:

1. Blood/plasma/VTM cell-free transfer medium: 200 μ l

2. Cell culture medium: at the rate of 1×10^7 cells

3. Homogenized tissue: add 30 to 50 mg of homogenized tissue in a volume of 400 microliters of tissue lysis buffer and homogenized completely. Finally add 200 microliters of the homogenized sample the next steps. Pay attention that incubation should not last more than 5 minutes at this stage.

Preparation of magnet nanoparticles: To prepare magnet nanoparticles, add 20 ml of RNase/DNase free distilled water included in the kit to the vial containing magnet nanoparticles and mix well.

RNA extraction steps

- Add 500 microliters of lysis buffer A to a RNase/DNase free microtube.
- Add the appropriate amount of the sample (for example, 200 microliters of blood or tissue lysate) to the microtube and mix vigorously for 15 seconds on the vortex and incubate at room temperature for 10 minutes. Vortex the samples every 5 minutes.
- Add 100 μ L of Lysis Buffer B to the microtube and mix vigorously for 5 seconds.
- Add 200 microliters of magnet nanoparticles to the microtube (be sure to shake the vial well before using the nanoparticles) and mix vigorously with a vortex and incubate for 5 minutes at room temperature.
- Place the microtube in a special rack and incubate for 30 seconds. The magnet particles that are attached to the RNA are moved towards the magnet in the rack. Therefore, by using the sampler, drain the liquid in the rack so that it does not come into contact with the magnet particles from inside the microtube.
- Remove the microtubes from the special rack and again add 300 microliters of lysis buffer A to the microtubes and vortex.
- Put the microtube in a special rack and use the sampler to drain the liquid in the rack so that it does not come into contact with the particles of the microtube.
- Remove the microtubes from the special rack and add 800 microliters of washing buffer to the microtubes and mix thoroughly by hand.
- Put the microtube in a special rack and use the sampler to drain the liquid in the rack so that it does not come into contact with the particles of the microtube.
- Add 100 microliters of RNase/DNase free distilled water, vortex and incubate for 2 minutes at room temperature. Make sure that the temperature of the distilled water is around 55 degrees Celsius for better RNA separation.
- Place the microtube in a special rack and incubate for 30 seconds. The liquid in the microtube contains pure RNA. Therefore, by using the sampler, transfer the existing liquid from the microtube to a sterile and RNase/DNase free microtube in such a way that it does not come into contact with the particles of the magnet and store it at -20°C .

<http://kpgene.ir>

09132926113

@karmaniaparsgene1

info@kpgene.ir

karmaniaparsgene

Troubleshooting:

- The current kit is highly quality controlled and the probability of encountering a problem is very low. However, in some cases, you may encounter problems using this kit. Refer to the following table to solve these issues:

Problem	Reasons	Solving
Low Yeild	Low cell count in the sample	✓ Use more samples
		✓ Add 50 μ L RNase/DNase free water
	Using an expired kit	✓ Although the components of the kits has more than 1 year stability, lysis buffer A and B had 6 months stability.
	Using the material contaminated with RNase	✓ Use the RNase/DNase free materials, You can prepare them from Karmania Pras Gene Company
	Using low amounts of lysis buffer	✓ Please follow the protocol
	Using expired magnets	✓ After adding the water, magnetic nanoparticles are active for 3 months
RNA is not appropriate for PCR	Presence of PCR inhibitor	✓ Do not change the washing step

Safety:

The lysing solution in the kit has a very alkaline and acidic Ph. On the other hand, the substances used in these solutions may in some cases be stimulating for the nervous system. Therefore, when working with these solutions, be sure to use gloves and in case of contact with the skin, wash with large amounts of water. In case of contact with the eyes, wash with large amounts of water and if there is no improvement, refer to medical centers.

KPG
Karmania Pars Gene