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Required items that are not included in the kit:

1. Sampler and sampler head
2. Centrifuge
3. ELISA reader containing a wavelength of 593 nm

Preparation of solutions:

- To prepare solution A in the 96-test kit, add 25 ml (12.5 ml in the 48-test kit) of distilled water to each vial. Note that due to the long-term instability of solution A after the addition of distilled water, two containers are provided in the kit to allow for two uses. Therefore, prepare solution A if needed.

Preparation of working solution:

- Mix the prepared solution A with acetic acid in a ratio of 1 to 3 and in the required amount. For example = 100 microliters of solution A with 300 microliters of acetic acid. (Solution A mixed with the prepared acid is stable for 72 hours in a dark place).
- Solution B and C should be mixed before testing. Accordingly, mix them in the ratio of 200/1. For example, if you need to test 4 samples, mix 995 µL of solution B with 5 µL of solution C. (The mixture is stable for only 48 hours and should be prepared immediately before testing). The solutions themselves are stable for 1 year if not mixed (the form provided in the kit).

Sample Preparation:

If using serum, plasma, or urine, no sample dilution is required.

For tissue samples: First, thoroughly homogenize 50 mg of tissue in 500 µL of sample dilution buffer (perform this step at a cold temperature; for example, pre-cool the dilution buffer by incubating it in the refrigerator for one hour). Then centrifuge at 5000 RPM for 5 minutes, and finally use the supernatant for the test.

Do not use the supernatant from cell culture samples, as it contains pH-sensitive colored reagents. Instead, first collect the cell sample from the culture medium, wash it with PBS, then fully lyse the cells in RIPA buffer specifically designed for oxidative stress tests (available from Karmania Pars Gene) using a sonicator. Centrifuge the cell lysate and use the supernatant for the test.

Note: If you do not have a sonicator, freeze-thaw the cell sample submerged in RIPA buffer 10 times, vigorously mixing each time.

Test Procedure:

1. Before starting the test, keep the kit components at room temperature for 20 minutes.
2. Add the ingredients according to the table below into 1.5 ml microtubes.

	Test	Negative Control	Standard	Blank
Sample	25 µL	25 µL	-	-
Distilled water	-	250 µL	25 µL	275 µL
Mix B and C solution	250 µL	-	250 µL	-

3. Incubate the tubes for 3 minutes at 37°C and then add 500 µl of the prepared solution A to all tubes except the blank.
4. Incubate all tubes for 10 minutes at 100°C and then centrifuge at 2500 rpm for 5 minutes.
5. Transfer 200 µL of the supernatant to the wells of the ELISA plate and evaluate its OD at 570 nm.
6. Use the following formula to calculate the catalase activity:

$$\text{Catalase Activity } \left(\frac{\text{KU}}{\text{L}} \right) = [S' - (S - M)] \times \frac{2.33}{T} \times \frac{V_t}{V_s}$$

S'=Absorbance of standard, S= Absorbance of Test,

M= Absorbance of negative control, VT=Total Volume and VS= Volume of sample

Given that the result of dividing 2.33 by 3 is 0.7676666, and also the result of dividing the final volume by the sample volume is 16, and the result of multiplying 0.7676666 by 16 is 12.282, so use the following formula:
Catalase activity (KU/L)= [S'-(S-M)]× 12.282

Note that the higher the catalase activity, the less color (green) will be produced because more solution C is consumed and does not react with solution A. This is why the standard, which does not contain catalase, reacts completely with solution A. Therefore, in cases where the sample color remains orange and the samples become uniformly orange, it indicates that solution C has been completely consumed and the catalase level in the sample is high. In these cases, dilute the sample by at least half. It is recommended to test one sample at different dilutions before performing all tests and then prepare the remaining samples at the appropriate dilution (OD: 0.5-1).

For easier data calculation, please download the Excel file designed by KPG Company using the barcode below:

Note that the OD of the numerical standard is about 0.6 to 0.8. If this number is lower, change the values for mixing solutions B and C accordingly. For example, if the absorbance obtained for the standard is 0.3, this ratio should be 1.50. Therefore, it is recommended to first perform a test on the standard and, if a suitable OD is obtained for the standard, perform the test on the samples.

Safety Instructions: The solutions used in this kit are hazardous to human tissue, hence, in cases of contact with skin, eyes, etc., wash with water and go to hospital for additional treatment.

Similar kits in the field of oxidants and antioxidants produced by Karmania Pars Gene:

1. Total Antioxidant Activity (FRAP) Testing Kit
2. Total Antioxidant Activity (DPPH) Testing Kit
3. Total Antioxidant Activity (ABTS) Testing Kit
4. Total Antioxidant Activity (CUPRAC) Testing Kit
5. MDA Activity Testing Kit
6. NO Level Testing Kit
7. Paraonase-1 Activity Testing Kit
8. Protein Carbonyl Level Testing Kit
9. SOD Activity Testing Kit
10. Acetylcholinesterase Activity
9. ...

Other products:

1. ELISA for measuring cytokines
2. DNA and RNA extraction by column, sedimentation and magnetic nanoparticles methods
3. Measurement of miRNAs
4. Products related to cell culture such as FBS, culture medium and cancer and mesenchymal stem cells
4. Microtubes Real-Time PCR devices
5. DNA and RNA purification columns

Contact us

+989132926113

WWW.KARMANIAPARSGEN.IR

karmaniaparsgene@gmail.com