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**Other materials needed:** Samplers, Centrifuge

**Sample:** This kit is designed to extract DNA from blood and tissue. Using this kit, you will be able to purify DNA from serum, plasma, VTM cell-free transfer medium, gram-negative bacteria and cell cultures, and homogenized tissue samples in addition to blood. For this, use the following volumes:

1. Blood/plasma/VTM cell-free transfer medium: 200  $\mu$ l
2. Cell culture medium: at the rate of  $1 \times 10^7$  cells
3. Homogenized tissue: add 30 to 50 mg of homogenized tissue in a volume of 400 microliters of homogenizing buffer (CN: KPG-KHBt). Then mix 200 microliters of the homogenized tissue with 400 microliters of lysis buffer A and after incubate for 5 minutes, continue the extraction steps. Pay attention that incubation should not last more than 5 minutes at this stage.

### Protocol:

- Add 500 microliters of lysis buffer B to the microtube provided in the kit.
- Add appropriate amounts of the sample (for example, 200 microliters of blood or all the tissue lysate) to the microtube and mix vigorously for 15 seconds on the vortex and incubate at room temperature for 10 minutes. Vortex the samples every 5 minutes.
- Add 100  $\mu$ L of Lysis Buffer C to the microtube and mix vigorously for 5 seconds.
- Add 200 microliters of magnet nanoparticles to the microtube (be sure to shake the vial well before removing the nanoparticles) and mix vigorously with a vortex and incubate for 5 minutes at room temperature.
- Place the microtube in a special magnetic rack and incubate for 30 seconds. The magnet particles attached to the DNA are moved towards the magnet in the rack. Therefore, using the sampler, drain the liquid.
- Remove the microtubes from the special rack and again add 300 microliters of lysis buffer B and 100 microliters of lysis buffer C to the microtubes and vortex.
- Put the microtube in a special rack and use the sampler to drain the liquid in the rack from inside the microtube.
- Remove the microtubes from the special rack and add 800 microliters of washing buffer to the microtubes and mix thoroughly by hand.
- Put the microtube in a special rack and use the sampler to drain the liquid in the rack.
- Add 70 microliters of DNase/RNase free distilled water, vortex and incubate for 2 minutes at room temperature. Make sure that the temperature of the distilled water is around 55°C for better DNA purification.
  - Place the microtube in a special rack and incubate for 30 seconds. The liquid in the microtube contains pure DNA. Therefore, using the sampler, transfer the existing liquid from the microtube to a sterile and DNase/RNase free microtube and store it at -20°C.

**Safety instructions:** The solutions used in this kit are dangerous for human tissue, so be sure to work with gloves and gloves, and in case of contact with skin, eyes, etc., wash with plenty of water and go to the hospital for additional treatment.

### Other DNA extraction kits and RNA produced by Karmania Pars Gene:

All the following DNA extraction kits are available in all formats, including column, sedimentary and magnet nanoparticles:

1. DNA extraction kits from blood and tissue by column
2. DNA extraction kits from gram-positive and negative bacteria
3. DNA extraction kit from mycobacterium
4. DNA extraction kit from Fungi
5. DNA extraction kit from virus
6. DNA extraction kit from HPV virus
7. DNA extraction kits from plant tissue

### RNA extraction kits:

All the following RNA extraction kits are available in both column and sedimentary methods:

1. RNA extraction kits from blood and tissue
2. RNA extraction kits from gram-positive and negative bacteria
3. RNA extraction kit from mycobacterium
4. RNA extraction kit from Fungi
5. The RNA extraction kit from the virus
6. RNA extraction kits from plant tissue

One of the most efficiently discharged kits:

### Plasmid extraction Kit

### Contact us

00989132926113

[karmaniaparsgene@gmail.com](mailto:karmaniaparsgene@gmail.com)

[WWW.KARMANIAPARSGEN.IR](http://WWW.KARMANIAPARSGEN.IR)

@karmaniaparsgene