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SOD Activity Assay Kit

Catalog Nos. KPG-SOD

The enzyme SOD (superoxide dismutase) is an antioxidant enzyme that plays an important role in protecting cells from free radical damage. This enzyme catalyzes the decomposition reaction of the superoxide radical (O₂⁻) and converts it into glass molecules and hydrogen peroxide, which are less toxic. Solution B in the Parsgene Carmania kit changes color in the buffered solution A by reacting with superoxide (O₂⁻) and is read at a wavelength of 420 nm. The SOD enzyme in the sample disrupts the reaction of solution B with superoxide (O₂⁻) by breaking down superoxide (O₂⁻), leading to a decrease in color. Therefore, the activity of this enzyme is directly related to the decrease in color.

➤ Kit components:

1. **Sample Reagent** (CN. KPG-PB7)
2. **SOD Solution A** (CN. KPG-SODA)
3. **SOD Solution B** (CN. KPG-SODB): 2 vials
4. **SOD Solution C** (CN. KPG-SODC)

Required items that are not included in the kit:

1. Sampler and sampler head
2. Centrifuge
3. ELISA reader containing a wavelength of 420 nm

Preparation of solutions:

- To prepare solution B in the 96-test kit, add 2.5 ml (1.25 ml in the 48-test kit) of solution C to solution B and mix well. Make sure that solution B is completely dissolved. Since solution B has a short shelf life after preparation, it is available in two separate containers in the kit so that the user can use the kit twice.
- The prepared solution B is stable for 72 hours in a dark place and in the refrigerator.

Sample Preparation: If you are using serum, plasma or urine, there is no need to dilute the sample. In cases where the sample is tissue, first thoroughly homogenize 50 mg of tissue in 500 µl of sample dilution buffer (this should be done at cold temperatures, for example, cool the dilution buffer by incubating it in the refrigerator for one hour). Then centrifuge at 5000 RPM for 5 minutes and finally use the supernatant to perform the test.

Test Procedure:

1. Before starting the test, keep the kit components at room temperature for 20 minutes.
2. Add 200 µL of solution A to the wells of the ELISA plate, according to the number of samples, and consider one well as a standard and add 210 µL of solution A to it.
3. Add 10 µl of the sample to the respective wells and after mixing immediately measure the absorbance of all wells at 420 nm. Save this absorbance as A1.
4. Add 50 µl of solution B to all wells and incubate for 10 minutes and then measure at 420 nm. Save this absorbance as A2.

Use the following formula to calculate SOD inhibitory activity:

$$\text{Inhibitory percent} = \left(1 - \frac{\text{Sample A2} - \text{Sample A1}}{\text{Standard A2} - \text{Standard A1}}\right) \times 100$$

Use the following formula to calculate enzyme activity units:

$$\text{SOD (U/mL)} = \left(\frac{\text{Inhibitory percent}}{50}\right) \times \text{Sample dilution}$$

Given that the sample dilution is 21 (10/210) and given that dividing 21 by 50 is 0.42, the following formula can be used:

$$\text{SOD Activity (U/ml)} = \text{Inhibitory percent} \times 0.42$$

For easier data calculation, please download the Excel file designed by KPG Company using the barcode below:

- Be careful if no color is observed in your sample or if it is very low, it is due to high enzyme activity in your sample and therefore you should dilute the sample with dilution buffer (for example, 10 microliters of sample with 20 microliters of dilution buffer). On the other hand, the absorbance (OD) in the standard should be between 0.8 and 1.3. Do not use cell culture supernatants directly due to the presence of pH-sensitive color indicators. Instead, cells must first be harvested from the culture medium and washed with PBS. Then, completely lyse the cells using a sonicator in RIPA buffer specialized for oxidative stress assays (available from Karmania Pars Gene). Centrifuge the cell lysate and use the resulting supernatant for testing.
- Note: If a sonicator is unavailable, subject the cell sample immersed in RIPA buffer to 10 cycles of freeze-thaw, ensuring vigorous mixing/vortexing during each cycle.
- Please note that for normalizing the data, you must also measure the total protein concentration (µg/mL) and divide the obtained SOD results by the total protein, reporting the data as U/µg. We recommend using the total protein assay kit from Karmania Pars Gene.

Safety Instructions: The solutions used in this kit are hazardous to human tissue, hence, in cases of contact with skin, eyes, etc., wash with water and go to hospital for additional treatment.

Similar kits in the field of oxidants and antioxidants produced by Karmania Pars Gene:

1. Total Antioxidant Activity (FRAP) Testing Kit
2. Total Antioxidant Activity (DPPH) Testing Kit
3. Total Antioxidant Activity (ABTS) Testing Kit
4. Total Antioxidant Activity (CUPRAC) Testing Kit
5. MDA Activity Testing Kit
6. NO Level Testing Kit
7. Paraoxonase-1 Activity Testing Kit
8. Protein Carbonyl Level Testing Kit
9. SOD Activity Testing Kit
10. Acetylcholinesterase Activity
9. ...

Other products:

1. ELISA for measuring cytokines
2. DNA and RNA extraction by column, sedimentation and magnetic nanoparticles methods
3. Measurement of miRNAs
4. Products related to cell culture such as FBS, culture medium and cancer and mesenchymal stem cells
4. Microtubes Real-Time PCR devices
5. DNA and RNA purification columns

Contact us
+989132926113

WWW.KARMANIAPARSGEN.IR
karmaniaparsgene@gmail.com