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Other materials needed: Samplers, Centrifuge, 1.5 mL tubes

Sample: This kit is designed to extract RNAs from blood and tissue. Using this kit, you will be able to purify RNAs from serum, plasma, VTM cell-free transfer medium, cell culture and homogenized tissue samples in addition to blood. For this, use the following volumes:

1. Blood/plasma/VTM cell-free transfer medium: 200 µl
2. Cell culture medium: at the rate of 1×10^7 cells
3. Thoroughly homogenize 30 to 50 mg of tissue in a volume of 400 µl of homogenizing buffer and finally use 200 µl of the sample for use in the next steps. Homogenization needs to be performed in 4 °C.

Preparation of solutions: All solutions are ready to use.

RNA extraction steps

- Add 500 microliters of lysis buffer A to a DNase/RNase free microtube.
- Add the appropriate amount of the sample (for example, 200 microliters of blood or 200 microliters of tissue lysate solution) to the microtube and mix vigorously for 15 seconds on the vortex and incubate at room temperature for 10 minutes. Vortex the samples every 5 minutes.
- Add 200 µL of Lysis Buffer B to the microtube and mix vigorously for 5 seconds.
- Add 200 µL of Lysis Buffer D to the microtube and mix vigorously for 5 seconds.
- Centrifuge the microtube at 12000-14000 RPM for 5 minutes and transfer the supernatant to another DNase/RNase free microtube.
- Add 300 microliters of cold precipitating buffer to the new microtube and mix gently by hand.
- Centrifuge the microtube at 12000-14000 RPM for 5 minutes and drain the supernatant completely.
- Add 500 microliters of washing buffer to the microtube and mix gently by hand.
- Centrifuge the microtube at 8000 RPM for 5 minutes and drain the supernatant completely and dry the remaining washing solution in the microtube by inverting it on a paper towel.
- Add 50 microliters of DNase/RNase free distilled water and incubate for 2 minutes at room temperature. Please note that the temperature of distilled water should be around 55 degrees Celsius to better isolate RNAs.

Safety instructions: The solutions used in this kit are dangerous for human tissue, so be sure to work with gloves and gloves, and in case of contact with skin, eyes, etc., wash with plenty of water and go to the hospital for additional treatment.

Other DNA extraction kits and RNA produced by Karmania Pars Gene:

All the following DNA extraction kits are available in all formats, including column, sedimentary and magnet nanoparticles:

1. DNA extraction kits from blood and tissue by column
2. DNA extraction kits from gram-positive and negative bacteria
3. DNA extraction kit from mycobacterium
4. DNA extraction kit from Fungi
5. DNA extraction kit from virus
6. DNA extraction kit from HPV virus
7. DNA extraction kits from plant tissue

RNA extraction kits:

All the following RNA extraction kits are available in both column and sedimentary methods:

1. RNA extraction kits from blood and tissue
2. RNA extraction kits from gram-positive and negative bacteria
3. RNA extraction kit from mycobacterium
4. RNA extraction kit from Fungi
5. The RNA extraction kit from the virus
6. RNA extraction kits from plant tissue

One of the most efficiently discharged kits: **Plasmid extraction Kit**

Contact us

00989132926113-00982191692296

karmaniaparsgene@gmail.com

WWW.KARMANIAPARSGEN.IR