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microRNA extraction kit from Blood and Tissue (CN: KPG-miRE)

A brief explanation of the microRNA extraction kit from the blood and tissues:

microRNAs are small molecules that are made of RNA and are able to bind to mRNAs and regulate the amount of translation. These molecules are involved in the pathogenesis of many diseases. Today, these molecules are very important in research and treatment. Column kits do not have the ability to purify small amounts of microRNAs. Therefore, most of the kits are designed in sedimentary form. Karmania Pars Gene microRNA extraction kit is able to purify microRNAs in sedimentary form with high capability and purity.

Kit component

Item	CN	volume
RNK Solution	KPG-DNK	50 mL (2 vial 25 mL)
Solution A	KPG-DNC	15 mL (2 vial 10 mL)
Precipitation solution	KPG-PS	15 mL
Carrier	KPG-Cr	1 mL
Homogenizing Buffer	KPG-Hbt	40 mL (2 vial 20 mL)
Washing buffer	KPG-WBp70	50 mL (2 vial 25 mL)
RNase/DNase free water	KPG-DW	5 mL
Microtube	KPG-M1/5	200 pcs

Material needs

-	Sampler tips
Centrifuge	Samplers

Sample: This kit is designed to extract microRNA from blood and tissue. Using this kit, you will be able to purify microRNA from serum, plasma, VTM cell-free transfer medium, and cell cultures, and homogenized tissue samples in addition to blood. For this, use the following volumes:

1. Blood/plasma/VTM cell-free transfer medium: 200 μ l
2. Cell culture medium: at the rate of 1×10^7 cells
3. Homogenized tissue: add 30 to 50 mg of homogenized tissue in a volume of 400 microliters of homogenizing buffer (CN: KPG-HBt) and finally use 200 microliters of the homogenized sample in next stage.

Protocol:

1. Add 500 microliters of RNK solution to a DNase/RNase free microtube.
2. Add the appropriate amount of the sample (for example, 200 microliters of blood, 200 microliters of tissue lysis solution) to the microtube and mix vigorously for 15 seconds on the vortex and incubate at room temperature for 10 minutes.
3. Add 200 μ L of Solution A to the microtube and mix vigorously for 5 seconds and incubate for five minutes.
4. Centrifuge the microtubes at 12.000 RPM for 5 minutes at 4°C.
5. Transfer the supernatant solution to a new microtube (be careful not to mix with the middle white layer) and add 300 microliters of cold sedimentation buffer and 10 microliters of Carrier to the microtube and mix gently by hand.
6. Centrifuge the microtubes at 12.000 RPM for 5 minutes at 4°C after a five-minute incubation.
7. Drain the supernatant completely by inverting the microtube and add 500 microliters of washing buffer and mix by hand and centrifuge at 8.000 RPM for 1 minute at 4°C.
8. Drain the supernatant completely by inverting the microtube and dehumidify it completely by placing it upside down on a paper towel.
9. Add 50 microliters of RNase/DNase free water and mix thoroughly. Incubate for 1 minute at room temperature. Please note that the temperature of RNase/DNase free water should be around 50 °C, so that microRNAs can be separated better.
10. Freeze the purified microRNAs in smaller microtubes (to prevent evaporation) at -20°C until use. Please note that the best mode of immediate use is for specific cDNA synthesis, but if needed for storage, microRNA can be stored for up to one month at -20°C and up to three months at -80°C.

Safety while using the kit

The solutions used in the kit have oxidizing and acidic properties. Avoid direct contact with skin and eyes. In case of contact with the mentioned tissues, wash with plenty of water and go to the nearest medical center.

Manufacture Address

Refsanjan, 20 km from Rafsanjan road to Kerman,
Rafsanjan special economic zone, food and drug area
09132926113
09135025983
Karmaniaparsgene@gmail.com
www.karmaniaparsgene.ir

Other DNA extraction kits and RNA produced by Karmania Pars Gene:

All the following DNA extraction kits are available in all three column, sedimentary and magnet nanoparticles:

1. DNA extraction kits from blood and tissue by column
2. DNA extraction kits from gram-positive and negative bacteria
3. DNA extraction kit from mycobacterium
4. DNA extraction kit from Fungi
5. DNA extraction kit from virus
6. DNA extraction kit from HPV virus
7. DNA extraction kits from plant tissue

RNA extraction kits:

All the following RNA extraction kits are available in both column and sedimentary methods:

1. RNA extraction kits from blood and tissue
2. RNA extraction kits from gram-positive and negative bacteria
3. RNA extraction kit from mycobacterium
4. RNA extraction kit from Fungi
5. The RNA extraction kit from the virus
6. RNA extraction kits from plant tissue

One of the most efficiently discharged kits: **Plasmid extraction Kit**