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Total Oxidant Status (TOS) Assay Kit

CN. KPG-TOS

Reactive oxygen and nitrogen species are produced during metabolic and physiological processes in the body. These compounds can cause harmful oxidative reactions that negatively affect organisms involved in antioxidant mechanisms (enzymatic and non-enzymatic). In certain circumstances, when the level of oxidants increases and the amount of antioxidants decreases, the body's defense system is unable to prevent oxidative damage, resulting in widespread disruptions in body function. Serum (or plasma) concentrations of different oxidant species can be measured separately, but these assay methods are usually time-consuming, expensive, and laborious, requiring complex laboratory techniques. Since it is not practical to evaluate individual oxidant molecules separately and their oxidative effects are cumulative, an index called total oxidant status (TOS) is used. This index measures the total amount of peroxides and oxidative activity present in the sample simultaneously and is used as a criterion for assessing oxidative stress. Due to the use of a single ready-to-use solution, this kit is very simple and easy to use compared to other kits on the market. The Total Oxidant Status Activity Measurement Kit produced by Karmania Pars Gene Company provides a standard, reliable, and simple method for assessing the total oxidant status in biological samples such as serum, plasma, urine, tissue, and cell lysates.

• Kit components:

1. Reagent TOS (KPG-RTOS)
2. TOS Dilution Buffer (KPG-PBSS2/7)
3. Standard Dilution buffer (KPG-PBSS6)
4. TOS Standard (KPG-TOSST)
5. ELISA Plate

Required items that are not included in the kit:

1. Sampler and sampler head
2. Centrifuge
3. ELISA reader containing a wavelength of 560 nm

Contact us

+989132926113

WWW.KARMANIAPARSGEN.IR

karmaniaparsgene@gmail.com

Preparation of solutions:

The total oxidant status reagent is in powder form and must be added to 50 ml of reagent dilution buffer in the 96-test kit (25 cc in the 48-test kit). The reagent is stable at room temperature for one hour. Therefore, prepare the reagent when needed and use the kit in less than one hour. If the reagent is stored in a fully oxidized state after testing, the color of the reagent will change over time and become purple. If needed again, first add concentrated HCL acid to the reagent until the reagent turns orange and then use again. Note that in this case, the sensitivity of the kit will decrease.

Sample preparation:

In cases where you use serum, plasma, or urine, there is no need to dilute the sample, but if the TOS level in the sample is high, the final color will turn gray instead of purple, so in these cases the sample must be diluted (usually one-half). In cases where the sample is tissue, first homogenize 50 mg of tissue in 500 μ L of standard dilution buffer (this should be done at cold temperature, for example, cool the dilution buffer by incubating it in the refrigerator for one hour) and then centrifuge at 10,000 RPM for 10 minutes. Finally, dilute the supernatant to one-fifth (e.g., 20 μ L of sample with 80 μ L of standard dilution buffer) with standard dilution buffer and then use in the test. Note that if the OD of the sample ultimately exceeds the standard of 2.5 mM/mL, the sample must be diluted further. Do not use the cell culture supernatant because it contains pH-sensitive color indicators. Instead, you must first collect the cell sample from the culture medium, wash it with PBS, then completely lyse the cells in RIPA buffer specifically for oxidative stress assays (you can obtain this from Karmania Pars Gene) using a sonicator. Centrifuge the cell lysate and use the supernatant for the assay. Please note that if you do not have a sonicator, you can subject the cell sample immersed in RIPA buffer to 10 cycles of freeze-thawing, with vigorous mixing each time.

Standard Preparation:

Add 5 cc of Standard Dilution Buffer to the vial containing the standard and mix well. This will prepare a 2.5 mM/mL standard solution. To prepare 0.83, 0.27, and 0.09 mM/mL solutions, follow the table below. The kit standard is stable for 1 month at room temperature.

Final concentration	Standard Dilution buffer	2.5 mM/mL standard
2.5 mm/mL	0 μ L	600 μ L
0.83 mm/mL	400 μ L	200 μ L
0.27 mm/mL	534 μ L	66 μ L
0.09 mm/mL	578 μ L	22 μ L

Test Procedure:

1. Before starting the test, keep the kit components at room temperature for 20 minutes.
2. Add 10 μ L of the sample and 10 μ L of each of the 2.5 to 0.09 mM standards to appropriate microtubes (e.g. 1.5 cc microtubes). Add 10 μ L of standard dilution buffer to one of the microtubes as a blank.
3. Add 400 μ L of Total Oxidation Status Reagent to all microtubes and mix well. Be careful to add the reagent to the sample microtube first, then to the standards and blank microtubes.
4. Mix the microtubes and incubate at room temperature for 4 minutes. Then transfer 200 μ L of the final solution to the ELISA plate and measure at a wavelength of 560 nm. Note that the color change from orange to purple indicates the reaction.

Calculations: Report data according to the standard curve and in millimolar per milliliter.

To calculate the total oxidative status in milligrams of protein, you need to measure the total protein in the sample. We recommend the Carmania Parsgene Total Protein Assay Kit because it is based on BCA technology and has much higher sensitivity than the Bradford method.

If calculating the total protein content, divide the data obtained by the total protein content and report the data as millimolar per microgram (mM/ μ g).

Safety Instructions: The solutions used in this kit are hazardous to human tissue, hence, in cases of contact with skin, eyes, etc., wash with water and go to hospital for additional treatment.

Similar kits in the field of oxidants and antioxidants produced by Karmania Pars Gene:

1. Total Antioxidant Activity (FRAP) Testing Kit
2. Total Antioxidant Activity (DPPH) Testing Kit
3. Total Antioxidant Activity (ABTS) Testing Kit
4. Total Antioxidant Activity (CUPRAC) Testing Kit
5. MDA Activity Testing Kit
6. NO Level Testing Kit
7. Paraonase-1 Activity Testing Kit
8. Protein Carbonyl Level Testing Kit
9. SOD Activity Testing Kit
10. Acetylcholinesterase Activity
9. ...

Other products:

1. ELISA for measuring cytokines
2. DNA and RNA extraction by column, sedimentation and magnetic nanoparticles methods
3. Measurement of miRNAs
4. Products related to cell culture such as FBS, culture medium and cancer and mesenchymal stem cells
4. Microtubes Real-Time PCR devices
5. DNA and RNA purification columns