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WWW.KARMANIAPARSGEN.IR

### Total protein Assay Kit Catalog Nos. KPG-ToP

The BCA Protein Quantitation Kit is a simple method for determining the concentration of total proteins. In this kit,  $\text{Cu}^{2+}$  binds to proteins and forms a color change, forming an intensely colored complex with bicinchoninic acid, which is read at a wavelength of 562 nm.

• **Kit Contents:**

Solution A (CN. KPG-TOPA)

Solution B (CN. KPG-ToPB)

Standard (CN. KPG-ToPS)

**Required items not included in the kit:**

Sampler and sampler head

Centrifuge

ELISA plate

ELISA reader containing a wavelength of 562 nm

**Preparation:**

- To prepare the **working solution**, mix solution A with solution B in a ratio of 1/50. For example, mix 10 microliters of solution A with 490 microliters of solution B.
- The prepared working solution is stable for 24 hours in a dark place.

**Preparation of standards:**

First prepare six sterile microtubes and label them as standards 1, 2, 3, 4, 5 and 6. Add 450  $\mu\text{l}$  of the buffer in which your protein is dissolved to the first microtube and 250  $\mu\text{l}$  to the remaining microtubes. Then add 50  $\mu\text{l}$  of the protein to the first microtube and mix well to prepare a 1000  $\mu\text{g/mL}$  standard. Then add 250  $\mu\text{l}$  of the first standard (1000  $\mu\text{g/mL}$ ) to the second standard and mix well to prepare a second standard with a concentration of 500  $\mu\text{g/mL}$ . Repeat this process until the final standard is prepared to prepare 250, 125, 5.62 and 25.31  $\mu\text{g/mL}$  concentrations.

**Test procedure:**

1. Before starting the test, keep the kit components at room temperature for 20 minutes.
  2. Add 25  $\mu\text{L}$  of sample and standards to separate microtubes and 25  $\mu\text{L}$  of distilled water as a control to a separate microtube.
  3. Add 225  $\mu\text{L}$  of working solution to all microtubes and incubate at 60°C for 60 minutes.
  4. Add 100  $\mu\text{L}$  of the solution to all wells and measure at 562 nm. Note that if you do not have a 562 nm wavelength, you can also read at 545 nm.
- To calculate the protein content of the sample, first subtract the OD of the blank from the OD of the samples and standards, and then calculate the protein content by drawing a standard graph.

**Safety Instructions:** The solutions used in this kit are hazardous to human tissue, hence, in cases of contact with skin, eyes, etc., wash with water and go to hospital for additional treatment.

**Similar kits in the field of oxidants and antioxidants produced by Karmania Pars Gene:**

1. Total Antioxidant Activity (FRAP) Testing Kit
2. Catalase Activity Testing Kit 11. MDA Activity Testing Kit
3. NO Level Testing Kit
4. Paraonase-1 Activity Testing Kit
5. Protein Carbonyl Level Testing Kit
6. SOD Activity Testing Kit
7. Acetylcholinesterase Activity

**Other products:**

1. ELISA for measuring cytokines
2. DNA and RNA extraction by column, sedimentation and magnetic nanoparticles methods
3. Measurement of miRNAs 10. Products related to cell culture such as FBS, culture medium and cancer and mesenchymal stem cells
4. Microtubes Real-Time PCR devices
5. DNA and RNA purification columns

**Contact us**  
**+989132926113**

[WWW.KARMANIAPARSGEN.IR](http://WWW.KARMANIAPARSGEN.IR)  
[karmaniaparsgene@gmail.com](mailto:karmaniaparsgene@gmail.com)