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DNA extraction kit from Saliva (CN: KPG-DNAsA)

A brief explanation of the DNA extraction kit from the Saliva:

To extraction of DNA from Saliva, use directly the sample. In the cases of tissue, you need to homogenize the tissue and, hence, in the cases of tissue you do not need solution A and Homogenizer buffer.

Kit component

Item	Catalogue number	volume
Lysis solution A	KPG-KLS	30 mL
Lysis solution B	KPG-KLSB	20 mL
Lysis solution C	KPG-KLSC	5 mL
Homogenizer Buffer	KPG-HBt	25 mL
Precipitation solution	KPG-PSES	13 mL
Washing buffer PI	KPG-WBpI	27 mL
Washing buffer PII	KPG-WBpII	27 mL
RNase/DNase free water	KPG-DW	5 mL
High absorbance column tube	KPG-Column	50 pcs
Collection Tube	KPG-CT	100 pcs

Material needs	
Distilled Water	Sampler tips
Centrifuge	Samplers

Sample: This kit is designed to extract DNA from Saliva. Using this kit, you will be able to purify DNA from blood, serum, plasma, VTM cell-free transfer medium, gram-negative bacteria and cell cultures, and homogenized tissue samples in addition to Saliva. For this, use the following volumes:

1. Saliva, Blood/plasma/VTM cell-free transfer medium: 200 μ L
2. Cell culture medium: at the rate of 1×10^7 cells
3. Homogenized tissue: add 30 to 50 mg of homogenized tissue in a volume of 400 microliters of homogenizing buffer (CN: KPG-HBt) and finally add 150 microliters of the homogenized sample in a new sterile tube and mix with 500 microliters of lysis solution A and after incubate for 5 minutes, continue the extraction steps. Pay attention that incubation should not last more than 5 minutes at this stage.

Protocol:

- Add 350 microliters of Lysis Buffer B to a DNase/RNase free microtube.
- Add an appropriate amount of the sample (for example, 200 microliters of Saliva, 200 microliters of 0.5 McFarland Gram-negative bacteria suspension, or the entire lysate solution obtained from the tissue or lysed bacteria) to the microtube and shake vigorously for 15 seconds. Mix on a vortex and incubate at 55 °C temperature for 20 minutes. Vortex the samples every 10 minutes.
- Add 50 μ L of Lysis Buffer C to the microtube and mix vigorously for 5 seconds.
- Add 250 microliters of cold precipitating buffer to the microtube and mix gently by hand. Be careful, you can mix Lysis Buffer C and precipitating buffer together and then add it to the microtube.
- Add 750 μ L of the contents of the microtube to the DNA extraction columns and centrifuge at 12000 RPM for 1 minute and 30 seconds. Make sure that according to the kit method, very good amounts of DNA are purified and you should not add extra sample volume to the column.
- Transfer the column to a new collection tube and add 530 μ L of wash buffer 1 and centrifuge at 8.000 RPM for 1 minute.
- Repeat the washing process with the washing buffer 2.
- Discard the collection tube content and centrifuge the column at 12.000 RPM for 1 minute.
- Transfer the column to a sterile DNase free 1/5 microtube and add 50 microliters of RNase/DNase free water and incubate for 1 minutes at room temperature. Make sure that the temperature of the distilled water is around 55°C for better DNA purification.
- Next, centrifuge the columns at 12/000 RPM for 2 minute. The resulting solution contains DNA. To increase the amount of extracted DNA, we suggest transferring the extracted DNA sample to the column and centrifuging again at 12.000 RPM for 1 minute.

Safety while using the kit

The solutions used in the kit have oxidizing and acidic properties. Avoid direct contact with skin and eyes. In case of contact with the mentioned tissues, wash with plenty of water and go to the nearest medical center.

Manufacture Address

Refsanjan, 20 km from Rafsanjan road to Kerman, Rafsanjan special economic zone, food and drug area
09132926113
09135025983
Karmaniaparsgene@gmail.com
www.karmaniaparsgen.ir

Other DNA extraction kits and RNA produced by Karmania Pars Gene:

All the following DNA extraction kits are available in all three column, sedimentary and magnet nanoparticles:

1. DNA extraction kits from blood and tissue by column
2. DNA extraction kits from gram-positive and negative bacteria
3. DNA extraction kit from mycobacterium
4. DNA extraction kit from Fungi
5. DNA extraction kit from virus
6. DNA extraction kit from HPV virus
7. DNA extraction kits from plant tissue

RNA extraction kits:

All the following RNA extraction kits are available in both column and sedimentary methods:

1. RNA extraction kits from blood and tissue
2. RNA extraction kits from gram-positive and negative bacteria
3. RNA extraction kit from mycobacterium
4. RNA extraction kit from Fungi
5. The RNA extraction kit from the virus
6. RNA extraction kits from plant tissue

One of the most efficiently discharged kits: **Plasmid extraction Kit**