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۰۹۱۳۲۹۲۶۱۱۳

WWW.KARMANIAPARSGEN.IR
karmaniaparsgene@gmail.com

Free Radical Scavenging Capacity Assay Kit based on DPPH CN. KPG-GPX96

The DPPH (2,2-diphenyl-1-picrylhydrazyl) method is one of the most widely used methods for determining the free radical neutralization capacity in plant and animal samples (such as serum, plasma, tissue, etc.) and also pharmaceuticals. Although this test can be performed on serum or plasma, milk, and some other samples that contain a large amount of protein or lipid with some changes in the working method and with some considerations. The DPPH method is briefly as follows: a sample of the desired substance is mixed with a DPPH solution. Then, using an ELISA reader, the decrease in DPPH concentration over time is measured. Using light absorption data, the antioxidant activity of the sample can be calculated. In such a way that antioxidant agents cause DPPH consumption, and the higher the antioxidant activity, the lower the color of the substrate.

Kit Contents:

1. Diluent Reagent A (KPG-DPTE8): 70 cc
2. Diluent Reagent B (KPG-DPSUBMETH): 15 cc
3. Substrate (KPG-DPSUBDP)
4. Proteinase K (KPG-PK): 1 cc
5. DPPH Standard (KPG-DPPHSTD)
6. ELISA Plate (96 Wells)

Required items not included in the kit: 1. Sampler and sampler head, 2. ELISA reader containing 517 nm wavelength.

Preparation of solutions: •

- To prepare the substrate, in 96 tests add 12.5 cc (6.5 cc in 48 tests kit) of Diluent Buffer B to the substrate container and mix well. Note that this solution should be stored in the dark and refrigerated and is stable for up to one week.
- To prepare the standard, add 64 cc of distilled water to the standard container. The prepared standard contains 50 μM (micromolar) of antioxidants. To prepare the first standard, 25 μM , mix 0.5 cc of the prepared 50 μM standard with 0.5 cc of Diluent Buffer B and continue this process to prepare 12.5 μM , 6.25 μM , and 3.12 μM .

Sample Preparation:

When using serum, plasma, or urine, there is no need to dilute the sample. However, if the protein concentration in the sample is high, the sample should be diluted (usually 1:2).

For tissue samples, first homogenize 20 mg of tissue in 500 μL of Dilution Buffer A thoroughly (this should be done at a cold temperature—for example, chill the dilution buffer by incubating it on ice for one hour). Then, centrifuge at 10,000 RPM for 10 minutes. Finally, dilute the supernatant 1:5 with Dilution Buffer A (e.g., 20 μL sample + 80 μL Dilution Buffer A) before use in the test.

Note: If the color reduction is not significant, the dilution factor should be decreased (e.g., 1:2). Therefore, before performing the test, always check one sample to determine the appropriate dilution. Do not use the cell culture supernatant for the test, as it contains pH-sensitive color indicators. Instead, first collect the cell sample from the culture medium, wash it with PBS, and then completely lyse the cells in RIPA buffer specifically intended for oxidative stress tests (you can obtain this from Karmania Pars Gene) using a sonicator. Centrifuge the cell lysate and use the resulting supernatant for the test. Note: If you do not have a sonicator, freeze and thaw the cell sample immersed in RIPA buffer 10 times, mixing vigorously each time.

Test Procedure:

Before starting the experiment, ensure all solutions are at room temperature.

Add 100 μL of the sample to a 1.5 mL microtube, followed by 100 μL of Dilution Buffer A and 10 μL of Proteinase K. Mix well and incubate at 55°C for 2 hours.

Transfer 100 μL of each sample and the standards (25, 12.5, 6.25, and 3.12 μmol) to ELISA wells. Designate one well as the blank and add 100 μL of Dilution Buffer B to it.

Add 50 μL of substrate to all wells. After mixing and incubating for 1 minute at room temperature, measure the absorbance of the samples, standards, and blank at 517 nm.

Subtract the blank absorbance from all readings and report the results in μmol based on the standard curve.

To express results as **percent inhibition**, use the following formula:

$$\text{Inhibition percent} = \frac{A_t}{A_b} \times 100$$

Where:

A_t = Absorbance of the sample at 517 nm

A_b = Absorbance of the blank at 517 nm

When evaluating a specific drug or compound, calculate the **IC₅₀** by selecting the dilution that yields 50% absorbance in this formula.

Notes:

If your sample does not contain protein, **Proteinase K addition and 2-hour incubation at 55°C are unnecessary**. Proceed directly to the next step after adding the sample and Dilution Buffer A.

- Since antioxidants decolorize the substrate, higher antioxidant levels result in less color.
- If the antioxidant concentration in your sample is high:
 - Either dilute the sample with Dilution Buffer A,
 - Or start the standard curve from **50 μM** and increase the substrate volume to **100 μL** .

Safety Instructions: The solutions used in this kit are hazardous to human tissue, hence, in cases of contact with skin, eyes, etc., wash with water and go to hospital for additional treatment.

Similar kits in the field of oxidants and antioxidants produced by Karmania Pars Gene:

1. Total Antioxidant Activity (FRAP) Testing Kit
2. Catalase Activity Testing Kit
3. MDA Activity Testing Kit
4. NO Level Testing Kit
5. Paraonase-1 Activity Testing Kit
6. Protein Carbonyl Level Testing Kit
7. SOD Activity Testing Kit
8. Acetylcholinesterase Activity
9. ...

Other products:

1. ELISA for measuring cytokines
2. DNA and RNA extraction by column, sedimentation and magnetic nanoparticles methods
3. Measurement of miRNAs 10. Products related to cell culture such as FBS, culture medium and cancer and mesenchymal stem cells
4. Microtubes Real-Time PCR devices
5. DNA and RNA purification columns

Contact us

+989132926113

WWW.KARMANIAPARSGEN.IR

karmaniaparsgene@gmail.com