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۰۹۱۳۲۹۲۶۱۱۳-۰۳۴۳۴۲۰۸۰۲۵

WWW.KARMANIAPARSGEN.IR

Ferric Reducing Ability of Plasma (FRAP) assay kit Catalog Nos. KPG-FRAP

Contents of the kit:

1. Solution A (CN. KPG-Sa)
2. Solution B (CN. KPG-Sb)
3. Solution C (CN. KPG-Sc)
4. Solution D (CN. KPG-Sd)
5. Solution E (CN: KPG-Se)
6. Standard (CN. KPG-StdFRAP)

Required items that are not included in the kit:

1. Sampler and sampler head
2. Centrifuge
3. ELISA reader containing a wavelength of 593 nm

Preparation of solutions:

- To prepare solution A, in a 96-test kit, add 1.25 ml (0.625 ml in a 48-test kit) of solution E to each bottle of solution A and mix well.
- To prepare solution B, in a 96-test kit, add 1.25 ml (0.625 ml in a 48-test kit) of deionized distilled water (DDW) to each bottle of solution B and mix well.
- To prepare solution C, in a 96-test kit, add 15 ml (7.5 ml in a 48-test kit) of solution D to each bottle of solution C and mix well.
- Note that solutions A, B and C are stable for one week after preparation. Therefore, in order for the user to be able to use the kit twice, there are 2 of these solutions in the kit. Therefore, prepare each one when needed.

Sample Preparation: In cases where you are using serum, plasma, or urine, sample dilution is not required. However, if the TAC activity level in the sample is very high, the final color will become too intense, and the OD will exceed 2.5, causing the samples to appear similar. Therefore, in such cases, the sample should be diluted. In cases where the sample is tissue, first thoroughly homogenize 50 mg of tissue in 500 μ L of dilution buffer (this should be done at a cold temperature; for example, pre-cool the dilution buffer by incubating it in the refrigerator for one hour). Then, centrifuge at 10,000 RPM for 10 minutes. Finally, dilute the supernatant one-fifth with the dilution buffer (e.g., 20 μ L of sample with 80 μ L of dilution buffer) and use it in the test. Note that if the color is not noticeable in the end, the dilution factor should be decreased (e.g., to 1/2). Do not use cell culture supernatant samples because they contain pH-sensitive-colored reagents. Instead, first collect the cell sample from the culture medium, wash it with PBS, then completely lyse the cells in RIPA buffer specifically for oxidative stress tests (you can obtain this from Karmania Pars Gene) using a sonicator. Centrifuge the cell lysate and use the supernatant for the test. Note: If you do not have a sonicator, freeze and thaw the cell sample immersed in RIPA buffer 10 times, mixing vigorously each time.

Preparation of working solution:

Solutions A, B and C should be mixed immediately before the test in a ratio of 1/10 to create a working solution. For example, if you need to examine 2 samples, prepare 300 μ L of working solution as follows: mix 25 μ L of solution A, 25 μ L of solution B and 250 μ L of solution C. The mixture is only stable for 30 minutes and should be prepared immediately before testing.

Preparation of standards:


Add 31 mL of deionized distilled water (DDW) to the standard bottle and mix well to obtain a 2 mM Fe^{2+} standard, and then make serial dilutions using the table below to create standards 1 to 6 using a 2 mM solution. create:

Standards are stable at room temperature for 3 months in the dark.

Std	Concentrations (Fe^{2+})	H_2O	Fe^{2+} 2mm
1	1 mM (1000 μ M)	500 μ L	500 μ L
2	0.7 mM (700 μ M)	650 μ L	350 μ L
3	0.5 mM (500 μ M)	750 μ L	250 μ L
4	0.3 mM (300 μ M)	850 μ L	150 μ L
5	0.2 mM (200 μ M)	900 μ L	100 μ L
6	0.1 mM (100 μ M)	950 μ L	50 μ L

Test method:

1. Keep the kit components at room temperature for 20 minutes before starting the test.
 2. Add 10 microliters of sample, standards and distilled water respectively to the wells of the ELISA plate corresponding to the sample, standards and blank.
 3. Add 140 μ L of working solution to all wells and mix for 30 seconds.
- Incubate the plate at 37°C for 5 minutes and then measure the OD of the samples along with the standards and the blank at a wavelength of 593 nm.

Sample standard curve. Please draw the standard curve yourself and do not use the curve.

Please note that for normalizing the data, you must also measure the total protein level, divide the TAC results by the total protein, and report the data as $\mu\text{M}/\mu\text{g}$. (For measuring total protein, we recommend using the Karmania Pars Gene total protein kit.)

Safety instructions:

The solutions used in this kit are hazardous to human tissue, therefore, in case of contact with skin, eyes, etc., wash with water and refer to hospital for additional treatment.

Contact us

00989132926113-00982191692296

karmaniaparsgene@gmail.com

WWW.KARMANIAPARSGEN.IR