Reverse Transcriptase (CN: KPG-RT)

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Components	Volume
RT buffer	4 µl
10mM dNTP mix	1 µl
Oligo-dT primers (10µM) Or Random Hexamer (10µM)	1 µl
Template RNA	5 µg
RT Enzyme	1 µl (1 unit/50 µl reaction)
RNase/DNase free distilled water up to	20 µl

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cDNA synthesis 45-65°C (47°C recommended)	60 minutes
Inactivation 95°C	10 minutes

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Reverse Transcriptase (CN: KPG-RT)

product introduction :

Many viruses use RNA as their genome. In the replication cycle, these viruses use an enzyme called Reverse Transcriptase (RT) to convert their RNA into DNA. In the laboratory, this enzyme is used to convert RNA to DNA, which is called Complementary DNA (cDNA) for many purposes. This enzyme lacks heat tolerance in the virus and in the wild model. On the other hand, in order to convert the purified RNA into cDNA, its coils must be opened under a certain temperature and exposed to the RT enzyme. Therefore, heat-resistant RT enzymes are very suitable for this purpose.

The wild form RT enzyme also has RNase property, which must be removed in the industrial form and used in the laboratory.

Application:


- Convert RNA to cDNA
- Used in cDNA synthesis kit
- Used in One-step mastermix

Storage buffer

Tris-HCl (pH 8.1) 20 mM, EDTA 0.1 mM, 2ME 8 mM, stabilizers and 58% glycerin

quality control:

1 microliter of genomic RNA was incubated with 10 units of RT in 50 microliter test buffer containing 1.5 mM MgCl₂ for 16 hours at 65°C, and the RNA was observed without degradation. Using the RT enzyme manufactured by Karmania Parsgen Company, cDNA synthesis was performed at temperatures of 45 to 65 degrees Celsius and the following results were obtained:

As it is clear in the picture, the synthesis of cDNA was carried out at temperatures of 45 to 65 degrees Celsius, and the enzyme tolerates temperatures above 60 degrees Celsius well, and the quality of conversion did not decrease.

Thermal stability:

10 units of the enzyme were incubated at 95 degrees Celsius for 30 minutes and then used in the PCR reaction, and at the end, the PCR product was analyzed by electrophoresis and confirmed the stability of the enzyme at high temperature.

Proposed protocol for cDNA synthesis

Use the following instructions as a general guideline:

Components

Volume

RT buffer	4 µl
10mM dNTP mix	1 µl
Oligo-dT primers (10µM) Or Random Hexamer (10µM)	1 µl
Template RNA	5 µg
RT Enzyme	1 µl (1 unit/50 µl reaction)
RNase/DNase free distilled water up to 20 µl	

Recommended temperature program:

cDNA synthesis	45-65°C (47°C recommended)
60 minutes	
Inactivation 95°C	10 minutes

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